

Versatile Powerful Mobile DIODE

Introducing

A.R.C. FOX 810 nm DIODE LASER

For Endo coagulation & Scleral Cyclophotocoagulation



Inexpensive
Versatile
Reliable
Flexible
Mobile
Powerful

Technical Data

Laser	Diode
Wavelengths - Output power*	810nm, > 5.5 Watt
Pulse length	1 ms to CW
Pulse interval	2 ms to 1 Sek
Aliming beam	640nm 2 mW
Dimensions	12 x 21 x 10 cm
weight	1.2 kg ready to use
Power requirement	Rechargeable batterie

Accessories

Endo Coagulation Probe	OD 19 gauge or OD 23 Gauge
Scleral Cyclophotocoagulation	Fiber 600 µ with foot plate
Protection Filter	Zeiss OPML, Leica
Slit Lamp adapter	Hang Streit Super View Hang Streit Clip-on Other Slit Lamps on request

* at Laserhead.

A.R.C. LASER



CARE  GROUP

DCR Surgery made Easy & Bloodless

Introducing 980 nm DCR LASER



Inexpensive
Versatile
Reliable
Flexible
Mobile
Powerful
9.1 w

Technical Data

Laser	Diode
Wavelengths - Output power*	980nm, > 9.1 Watt
Pulse length	1 ms to CW
Pulse interval	2 ms to 1 Sec
Aiming beam	640nm 2 mW
Dimensions	12 x 21 x 10 cm
Weight	1.2 kg ready to use
Power requirement	Rechargeable batteries

* at Laserhead.

A.R.C. LASER



Intra Ocular Care Pvt.Ltd.

B, South-West Compound, Near ITC Green, Kolkata - 700 016, India, 700016.
Ph: 91-361-2651230/2651231, F: 91-361-2651232, MO: 9830011411, P: 1-26-1230/14
Email: intraocularcare@iocare.com, iocare@iocare.com, www.intraocularcare.com, www.iocare.com